



O-TEM250 with Trolley

Description

These Open TEM-cells are well suited for immunity testing of small objects according to European (CE) and automotive standards (SAE J1113-25) or for biological experiments. The advantage of these O-TEM cells is that they are open and it is very easy to control the functions of the EUT (Equipment Under Test). The applications are for instance the immunity testing of watches, pagers, telephones or PCB's.

In comparison with other closed TEM-cells, the price is low. The field decreases rapidly outside the Open TEM-cells (approx. 33 dB at 1 meter) and it is therefore possible to use an Open TEM-cell in ordinary facilities. Another interesting application is the calibration of field probes as the generated field inside the Open TEM-cell is very exact. Open TEM-cells are the most precise structures for field calibrations.

Technical specifications	O-TEM 250	O-TEM 500	O-TEM 1000 *	O-TEM 3000 *
Frequency range	DC - 250 MHz	DC - 500 MHz	DC - 1 GHz	DC - 3 GHz
Height under plate	30.0 cm	14.7 cm	7.3 cm	2.5 cm
Max. input power	1.5 kW	1 kW	750 W	400 W
Maximum field	800 V/m	1.5 kV/m	2.6 kV/m	5.6 kV/m
Field for a 25 W amplifier	105 V/m	215 V/m	480 V/m	1400 V/m
Dimensions (LxWxH)	250 x 114 x 65 cm	106 x 51 x 32 cm	54 x 45 x 16.8 cm	44 x 18 x 8 cm
Weight	50 kg	18 kg	3.5 kg	3 kg
Field precision	± 5 %			
Connector	type N 50 Ω, Optional : 7/16" or EIA 7/8" 50 Ω			
Cell impedance	50 Ω			
Wave impedance	377 Ω			
VSWR	< 1.2 (<1.9 for the O-TEM3000)			
Options	Dummy Load, Signal and power supply filters,			
	* Can be equipped with a test setup for IC or PCB acc. to SAE J1752-3			

Specifications subject to change without notice.

Other models are available on request

Sales Partner:



ABSOLUTE EMC LLC.
Covering sales in North America
United States, Mexico, & Canada

absolute-emc.com
Phone: 703-774-7505
info@absolute-emc.com